

# The SAFER*plus*

## The smallest & safest pass-through AER on the market

The SAFER*plus* is fully compatible with Sterilox - the only disinfectant on the market with no COSHH implications. Sterilox is safe for both patients and users, with no health and safety risks.

### Features

- **Quick**
- **Small**
  
- Fully compatible with Sterilox high level disinfectant and a wide range of Sterilox endoscopy products
- Pass-through design
- Can be used single sided
- Modular design – available in 1 or 2 endoscope design
- Dual display
- Available with or without Reverse Osmosis
- Fully meets your COSHH obligations
- Protects staff & patients from exposure to hazardous chemicals
- Reliable – built on proven technology
- Full tracking & traceability
- Small footprint
- Designed to meet EN15883-4 and HTM2030



**sterilox**  
endoscopy  
a PuriCore solution

# The SAFER<sub>plus</sub>

## Technical Features - The SAFER<sub>plus</sub>

Technical Features	
Disinfectant	Sterilox (HOCl)
Scopes per cycle	1 or 2 scope – re-processed in approximately 23-27 mins 5 min disinfectant contact time!
Rinse Water	Sterilox Rinse Water or Reverse Osmosis
Chamber Volume	35 litres
Documentation	Hard print and data capture – Integral Printer
Display	Dual display on clean & dirty side
Scope Compatibility	All flexible endoscopes
Overall Size	2027mm(H) x 1062mm(W) x 841mm(D)
Construction	Stainless Steel
Power Requirements	230vAC+10% 5A
Other Specifications	Designed by HTM2030 accredited engineers.

For further information please contact a member of the Sterilox team at the number below or enter your details on our website [www.puricore.com](http://www.puricore.com)



Distributed by; Sterilox endoscopy

PuriCore International Ltd  
Wolseley House, Staffordshire Technology Park,  
Beaconside, Stafford ST18 0GA. UK

Tel: +44 (0) 1785 782420

Fax: +44 (0) 1785 782421

Or visit us at: [www.puricore.com](http://www.puricore.com)

Part No. 0001-0428 Revision – A.

**PuriCore**  
Pure Science. Pure Life.